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#### Introduction

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**Sonicare For Kids** rechargeable sonic toothbrush

**PHILIPS** 


### www.philips.com/Sonicare


©2015 Koninklijke Philips N.V. (KPNV). All rights reserved Sonicare and the Sonicare logo are trademark of Philips Oral Healthcare Inc. and/or KPNV

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#### **ENGLISH**

#### Introduction

Congratulations on your purchase and welcome to Philips! To fully benefit from the support that Philips offers, register your product at www.philips.com/welcome

You and your child now own a Sonicare connected toothbrush. The Bluetooth connected toothbrush interacts with a fun brushing mobile device application to keep your child engaged while they learn proper brushing techniques.

### Sonicare For Kids (Fig. 1)

- 1 Hygienic travel cap
- 2 Brush head
- 3 Handle with soft grip
- 4 Power on/off button with battery charge indicator
- 5 Removable panel with replaceable stickers
- 6 Brushing mode button with mode lights
- 7 Charger
- 8 Panel 9 Sticker

# **IMPORTANT SAFEGUARDS**

### READ ALL INSTRUCTIONS BEFORE USE

#### **DANGERS**

- To reduce the risk of electrocution:
- Do not place or store the product where it can fall or be pulled into a tub or sink, Do not place the product in or drop it into water or any other liquid.
- Do not reach for a product that has fallen into water. Unplug immediately
- 4 Do not use the product while bathing.

#### WARNINGS

- To reduce the risk of burns, electrocution, fire or injury to persons: Discontinue use if this product appears damaged in any way (brush head, handle, charger). This product contains no serviceable parts. Refer to chapter 'Warranty and support' if your Sonicare no longer works properly or needs repair.
- Never use the charger if it has a damaged cord or plug.
- Keep the cord away from heated surfaces. Do not use the charger outdoors.
- Do not use the charger if dropped into water.
- This product is designed to clean your teeth, gums and tongue only. Use this product only for its intended use as described in this booklet. Discontinue use of this product and contact a physician/dentist if
- discomfort or pain is experienced. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for
- their safety. 8 Young children learning to brush should be supervised to ensure that they do not play with the appliance.
- Do not use attachments other than those recommended by the
- 10 If your toothpaste contains peroxide, baking soda or other bicarbonate (common in whitening toothpastes), thoroughly clean the brush head and the handle with soap and water after each use. Otherwise the plastic may crack
- 11 Follow carefully all the safety and warning instructions provided by the manufacturer when you use other devices near water and wet locations. Use the instructions provided by the manufacturer to address any food or liquid spills.
- 12 Do not charge your phone or tablet in the bathroom.

### **MEDICAL WARNINGS**

- Consult your dental professional before you use this product if you have had oral or gum surgery in the previous 2 months.
- Contact your dental professional if excessive bleeding occurs after using this product or bleeding continues to occur after 1 week of use.
- If you have a pacemaker or other implanted device contact your physician or the device manufacturer prior to use.
- This Philips appliance complies with all applicable standards and regulations regarding exposure to electromagnetic fields.
- Consult your physician prior to using the Sonicare if you have medical

## **NOT A TOY**

CHECK BRUSH HEADS REGULARLY, CRACKS FROM DROPPING CAN LEAD TO A CHOKING HAZARD.

# **SAVETHESE INSTRUCTIONS**

# Applying the replaceable panel

The Sonicare For Kids comes with 1 replaceable panel.

1 Remove the adhesive backing from the panel.

2 Apply the panel to the white section on the front of the handle (Fig. 2). Note: Make sure you line up the cut-out areas on the panel with the buttons on the front of the handle.

3 Apply your favorite sticker to the panel as shown on the sticker sheet (Fig. 3).

Reapplying a panel that has been removed

- 1 Rinse the panel in warm water; you may also use a mild detergent to
- 2 Let the panel air dry. 3 Reapply the panel to the white section on the front of the handle.
- Attaching the brush head

1 Align the brush head so the bristles point in the same direction as

- the front of the handle (Fig. 4). 2 Firmly press the brush head onto the metal shaft until it stops.
- Note:There is a small gap between the brush head and the handle.

# **Charging your Sonicare**

- 1 Plug the charger into a live electrical outlet.
- 2 Place the handle on the charger (Fig. 5).
- ▶ The battery charge indicator in the power on/off button flashes and the handle beeps twice to indicate that the toothbrush is charging. When the Sonicare is fully charged, the battery charge indicator stops flashing and remains on.

Note: If the battery charge of your Sonicare is low, you hear beeps and the battery charge indicator flashes rapidly for 30 seconds after the brushing cycle. Note:To keep the battery fully charged at all times, you may keep your Sonicare

on the charger when not in use. When the battery is fully charged, the charger draws minimal energy from the electrical outlet. It takes at least 24 hours to fully charge the battery.

#### ■ Using the Sonicare For Kids Bluetooth connectivity

Before you start the connectivity procedure, make sure your mobile device runs on:

- Apple iOS 7 or higher. Bluetooth 4.0 or higher.
- Note: For the latest information on device compatibility, please visit the App

4235\_020\_9982\_1\_A6\_Leaflet\_v1.indd 2

Store or Google Play.

Connecting your toothbrush

- 1 Download the Philips Sonicare For Kids app from the App Store or
- 2 Confirm your smartphone or tablet Bluetooth is on before you open the app. Check the settings on your smartphone or tablet to turn on
- 3 To connect the handle to the app, turn on the handle by pressing either the power on/off button or the brushing mode button. Make sure the handle is not on the charger when you turn it on.

- 4 Open the Philips Sonicare For Kids app.
- 5 Follow the app instructions displayed on your smartphone or tablet. The app automatically connects with the toothbrush. If the toothbrush does not connect, follow the troubleshooting steps provided in the app.
- 6 Keep the smartphone or tablet close to your toothbrush handle when you use your toothbrush. Make sure the smartphone or tablet is safe in a dry place.

### **Brushing instructions**

- 1 Wet the bristles and apply a small amount of toothpaste (about the size of a pea).
- 2 Place the toothbrush bristles against the teeth at a slight angle (45 degrees), where the teeth and the gumline meet. This is to ensure that the whole tooth surface is brushed (Fig. 6).
- 3 Press the power on/off button to turn on the Sonicare.
- 4 Apply light pressure to maximize Sonicare's effectiveness and let the Sonicare toothbrush do the brushing for you (Fig. 7).
- Maintain light pressure and gently brush the outer and inner surfaces of all top and bottom teeth while you keep the toothbrush at an angle such that the toothbrush bristles touch the area where the teeth and the gumline meet.
- 6 The inner and outer surfaces of the upper and lower back teeth are hard to reach. Spend time brushing these areas, apply light pressure and maintain the angle of the toothbrush bristles to the gum line
- 7 Brush the chewing surfaces of all the upper and lower teeth.
- 8 While you maintain the angle of the bristles to the gum line, gently move the brush head slowly across the teeth in a small back and forth motion so the longer bristles reach between the teeth. Continue this motion throughout the brushing cycle (Fig. 9).

Note:To make sure you brush evenly throughout the mouth, divide the mouth into 4 sections using the KidPacer feature (see chapter "Features").

- 9 Begin brushing in section 1 (outside top teeth) a series of musical notes indicate it is time to move to section 2 (inside top teeth). At the next series of notes move to section 3 (outside bottom teeth). The next series of notes indicate it is time to move to section 4 (inside bottom teeth). The toothbrush will automatically turn off at end of the brushing cycle (Fig. 10).
- Your Sonicare is safe to use on: Braces (brush heads wear out sooner when used on braces) Dental restorations (fillings, crowns, veneers)

### **Brushing modes**

The toothbrush automatically starts in the default high mode.

1 Press the brushing mode button to toggle between modes (Fig. 11). The green mode light indicates the selected mode.

High mode: Full-strength mode designed for older/independent brushers. Low mode: Reduced power mode used for training and introducing younger children to the sonic experience.

#### Features **KidTimer**

- The KidTimer plays a congratulatory sound sequence at the end of the pre-determined brushing time. Additionally, the KidTimer automatically turns off the handle when the brushing cycle is complete. When in high mode, the KidTimer is set at 2 minutes.
- When in low mode, the KidTimer is set at 1 minute and slowly increases over time. The progression of time helps young children reach the dentist recommended 2 minutes of brushing.

#### Note: Instruct your child to brush until they hear the congratulatory tune and the toothbrush turns off.

### **KidPacer**

The KidPacer uses a short series of tones to alert the brusher to move to the next brushing quadrant of the mouth. This ensures that the child's teeth gets a thorough, overall clean (see "Brushing instructions").

14 brushings.

- The Sonicare For Kids comes with the Easy-start feature activated. The Easy-start feature gently increases the power over the first

Note: Each of the first 14 brushings must be at least 1 minute in length to properly advance through the Easy-start ramp-up cycle. Deactivating or activating the Easy-start feature

1 Attach the brush head to the handle

2 Place the handle in the plugged-in charger. To deactivate Fasy-start: Press and hold the power on/off button for 2 seconds. You hear 1 beep to indicate that the Easy-start feature has been deactivated.

To activate Easy-start: Press and hold the power on/off button for 2 seconds. You hear 2 beeps to indicate that the Easy-start feature has been activated.

Note: Using the Easy-start feature beyond the initial ramp-up period is not recommended and reduces Sonicare's effectiveness in removing plaque.

### Do not clean brush heads, handle, charger, charger cover or replaceable panel in the dishwasher.

Toothbrush handle 1 Remove the brush head and rinse the metal shaft area with warm

water (Fig. 12). Do not push on the rubber seal on the metal shaft with sharp objects,

as this may cause damage.

2 Use a damp cloth to wipe the entire surface of the handle.

- 1 Rinse the brush head and bristles after each use (Fig. 13).
- 2 Remove the brush head from the handle and rinse the brush head connection at least once a week with warm water.

- 1 Unplug the charger before cleaning.
- 2 Use a damp cloth to wipe the surface of the charger.

- If you are not going to use the Sonicare for an extended period of time, unplug the charger from the electrical outlet, clean it and store it

# in a cool and dry place away from direct sunlight.

Replacement

# **Brush head**

- Replace Sonicare brush heads every 3 months to achieve optimal results. Use only Sonicare For Kids replacement brush heads.
- Disposal This appliance contains a rechargeable battery which must be disposed
- of properly (Fig. 14). Contact your local town or city officials for battery disposal information. You can also call 1-800-8-BATTERY or visit **www.rbrc.com** for
- battery drop-off locations. For assistance, visit our website **www.philips.com/support** or call

#### 1-800-243-3050 toll free. Removing the rechargeable battery

cap and the handle (Fig. 15).

### Only remove the rechargeable battery when you discard the appliance. Make sure the battery is completely empty when you remove it.

- To remove the rechargeable battery, you need a flat-head (standard) screwdriver. Observe basic safety precautions when you follow the procedure outlined below. Be sure to protect your eyes, hands, fingers, and the surface on which you work.
- 1 To deplete the rechargeable battery of any charge, remove the handle from the charger, turn on the Sonicare and let it run until it stops. Repeat this step until you can no longer turn on the Sonicare. 2 Insert a screwdriver into the slot located in the bottom of the handle

and turn counterclockwise until you see a gap between the bottom

Insert the screwdriver into the gap and pry the bottom cap out of the handle (Fig. 16). 4 Hold the handle upside down and push down on the shaft to release

the internal components of the handle (Fig. 17).

5 Insert the screwdriver under the circuit board, next to the battery connections, and twist to break the connections. Remove the circuit board and pry the battery from the plastic carrier (Fig. 18). The rechargeable battery can now be recycled or disposed of and the rest of the product discarded appropriately.

## **■** Warranty and support

Philips guarantees its products for two years after the date of purchase. Register your product at **www.philips.com/support**. Defects due to faulty materials and workmanship will be repaired or replaced at Philips' expense provided that convincing proof of purchase in the qualifying period is provided. Philips recommends use of genuine Sonicare brush heads for optimal performance. Use of other brush heads may damage your Sonicare toothbrush and limit your warranty protection. Contact our Consumer Care Center at 1-800-682-7664 (North America), outside North America contact your local Philips Consumer Care Center. Internet information: www.sonicare.com (North America) or

#### www.philips.com/support (outside North America). **WARRANTY EXCLUSIONS**

What is not covered under warranty:

- Brush heads. Replaceable panels.
- Damage caused by misuse, abuse, neglect, alterations or unauthorized
- Normal wear and tear, including chips, scratches, abrasions, discoloration or fading.
- Replaceable stickers.

#### **IMPLIED WARRANTIES**

ANY IMPLIED WARRANTIES, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED TO THE DURATION OF THE EXPRESS WARRAN TIES SET FORTH ABOVE, IN SOME STATES LIMITATIONS ON DURA-TION OF IMPLIED WARRANTIES DO NOT APPLY.

#### **LIMITATION OF REMEDIES**

IN NO EVENT SHALL PHILIPS OR ANY OF ITS AFFILIATED OR SUBSIDIARY COMPANIES BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED UPON BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, TORT, OR ANY OTHER LEGAL THEORY, SUCH DAMAGES INCLUDE, WITHOUT LIMITATION, LOSS OF SAVINGS OR REVENUE; LOSS OF PROFIT; LOSS OF USE; THE CLAIMS OF THIRD PARTIES INCLUDING, WITHOUT LIMITATION, DENTISTS AND DENTAL HYGIENISTS; AND COST OF ANY SUBSTITUTE EQUIPMENT OR SERVICES. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES.

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